



Gulden Ophthalmics
225 Cadwalader Avenue
Elkins Park, PA 19027-2020, USA
Web: www.guldenophthalmics.com
Email: info@guldenophthalmics.com
Phone: +1(215) 884-8105



Qarad B.V.
Flight Forum 40
5657 DB Eindhoven
The Netherlands



The filters are marked according to the left column (“Density”). The relationship between density and standard log units is as follows:

<u>Density</u>	<u>Log Units</u>
03	0.3
06	0.6
09	0.9
12	1.2
15	1.5
18	1.8

Blotch marks may appear due to the interference phenomena when gel filters in proximity to the plastic sheet cause an air pocket. These do not affect the light transmission of the filter and therefore the function is unaffected. The air pockets may come and go with any change of temperature and humidity.