

Ellis Magnetic Foreign Body Remover

Instructions For Use:

Purpose:

Magnetic Foreign Body Removers are utilized to remove foreign bodies from a patient's eye. The magnetization allows metal pieces to be removed while reducing the rates of irritation and corneal tearing.

Cleaning:

The Ellis Magnetic Foreign Body Remover should be disinfected with a 5-10 minute soak in one of following disinfectants: 3% hydrogen peroxide, or 70% isopropyl alcohol prior to use.

Warning:

Product creates small magnetic field.

Do not place into contact with electronic devices.



225 Cadwalader Avenue, Elkins Park, PA 19027
(800) 659-2250 Fax (215) 884-0418 e-mail: info@guldenophthalmics.com
www.guldenophthalmics.com

Rev. EDA 1108 || 6/18/2021

Dix Magnetic Foreign Body Remover

Instructions For Use:

Purpose:

Magnetic Foreign Body Removers are utilized to remove foreign bodies from a patient's eye. The magnetization allows metal pieces to be removed while reducing the rates of irritation and corneal tearing.

Cleaning:

The Dix Magnetic Foreign Body Remover should be disinfected with a 5-10 minute soak in one of following disinfectants: 3% hydrogen peroxide, or 70% isopropyl alcohol prior to use.

Warning:

Product creates small magnetic field.

Do not place into contact with electronic devices.



225 Cadwalader Avenue, Elkins Park, PA 19027
(800) 659-2250 Fax (215) 884-0418 e-mail: info@guldenophthalmics.com
www.guldenophthalmics.com

Rev. EDA 1108 || 6/18/2021